

MR100 Plus™

MANUAL RESUSCITATOR



Silicone mask

Star lumen oxygen tubing

Patient valve

- 40/60 cmH₂O safety valve (for infant / adult) releases excessive pressure at fixed level
- 360° rotatable design to adjust mask for any direction



Compressed bag

- Rapid recoil for optimum ventilation
- Silicone bag with high transparency
- Bag volume 1,500 ml, 550 ml, 280 ml



All-in-1 intake assembly

- Reduced space and weight design for supplying oxygen efficiently
- Lighter and thinner reservoir bag for lowering flow resistance



Reservoir bag

- O₂ reservoir bag volume: 2,000 ml, 1000 ml, 600 ml

Housing	Compressed Bag	Reservoir Bag		Mask		Oxygen Tube	Valve
PSF	Silicone	Bag: Silicone	Coupling: PSF	Cushion: Silicone	Housing: PSF	PVC	Silicone

Order information

REF	Mask Size	Pop-Off Valve (cmH ₂ O)	Bag Capacity (ml)	Reservoir Bag (ml)	Oxygen Tube (cm)	QTY (boxes/carton)
AR0056	#4	60	1500	2000	200	12
AR0064	#5	60	1500	1000	N/A	12
AR0060	N/A	60	1500	N/A	N/A	12
AR0044	#3,4,5	N/A	1500	2000	200	12
AR0059	#5	N/A	1500	2000	200	12
AR0063	N/A	N/A	1500	N/A	N/A	12
AR0045	#1,2,3	40	550	1000	200	12
AR0057	#3	40	550	1000	200	12
AR0061	N/A	40	550	N/A	N/A	12
AR0046	#0,1	40	280	500	200	12
AR0058	#0	40	280	500	200	12
AR0062	N/A	40	280	600(PVC)	N/A	12



Accessories

Image	REF	Description	QTY
	AM0004	Sil.Flex Mask, #0	50 pcs/carton
	AM0005	Sil.Flex Mask, #1	
	AM0016	Sil.Flex Mask, #2	
	AM0023	Sil.Flex Mask, #3	
	AM0024	Sil.Flex Mask, #4	
	AM0086	Sil.Flex Mask, #5	
	AM0001	Guedel Airway 40 mm, Pink	200 pcs/carton
	AM0046	Guedel Airway 50 mm, Blue	
	AM0047	Guedel Airway 60 mm, Black	
	AM0048	Guedel Airway 70 mm, White	
	AM0005	Guedel Airway 80 mm, Green	
	AM0006	Guedel Airway 90 mm, Yellow	
	AM0049	Guedel Airway 100 mm, Red	
	AM0002	Guedel Airway 110 mm, Orange	
	AM0003	Guedel Airway 120 mm, Purple	
	AM0074	Bermen Airway 40 mm	200 pcs/carton
	AM0050	Bermen Airway 50 mm	
	AM0051	Bermen Airway 60 mm	
	AM0052	Bermen Airway 70 mm	
	AM0053	Bermen Airway 80 mm	
	AM0007	Bermen Airway 90 mm	
	AM0054	Bermen Airway 100 mm	
		AR0003	
AR0005		Disposable Diverter, 30/26 mm	
	AF0009	Filter, 30F, 30M	50 pkg/box
	AM0008	Filter, 30F, 22/15 mm	

Image	REF	Pressure range	Connector size	Material
<p>Reusable Adujustable PEEP Valve</p>	AV0003	5-20 (cmH ₂ O)	30F, 22M, 22/15 mm	PSF
	AV0016	2-10 (cmH ₂ O)	22/15 mm	PSF
	AV0017	5-20 (cmH ₂ O)	22/15 mm	PSF
	AV0018	5-20 (cmH ₂ O)	22M	PSF
	AV0035	2-10 (cmH ₂ O)	30F	PSF
	AV0037	5-20 (cmH ₂ O)	30F	PSF
<p>Disposable Adujustable PEEP Valve</p>	AV0002	5-20 (cmH ₂ O)	22/19 mm	PC
	AV0004	5-20 (cmH ₂ O)	30F	PSF
	AV0019	2-10 (cmH ₂ O)	30F	PC
	AV0092	5-20 (cmH ₂ O)	30F	PC

AP2 UG (haftungsbeschränkt), Erich-Ollenhauer-Str.12, 51377 Leverkusen
 Handelsregister: HRB 111772 Registergericht: Amtsgericht Köln
 Vertreten durch: Paul Stamirski, Anton Zirnow
 Telefon: +49 (0) 214 690 58 414, Email: info@ap2.tech